

Intercomparisons of the Harvard Lyman alpha hygrometer and ICOS isotopic water instrument with the CFH and MLS instruments: Implications of recent results

Elliot M. Weinstock, J. B. Smith, R. Lockwood, D. S. Sayres, T. F. Hanisco, E. J. Moyer, J. M. St. Clair, and J. G. Anderson, Department of Chemistry and Chemical Biology, Harvard University, Cambridge MA 02138


W. G. Read, Jet Propulsion Laboratory, California Institute of Technology, Pasadena, CA

With acknowledgements and thanks to Holger Voemel, Bob Herman, Chris Webster for the use of their data.


Questions to be explored



- Are intercomparison data from CRAVE and AVE-WIIF consistent?
- What have we learned from CRAVE regarding the accuracy of in situ water instruments needed for Aura satellite validation, especially regarding the previously observed systematic differences between the frost point hygrometer and in situ aircraft instruments?
- How do MLS version 1.5 and version 2 compare with in situ water vapor measurements?

AVE-WIIF 20050703



AVE-WIIF 20050705





Water measurements during AVE-WIIF

Harvard Lyman- α
water vapor


Pallet 3
ICOS water vapor isotopes

Pallet 2
HOxotope
water vapor
isotopes

JLH multipass laser IR
water vapor

Pallet 4
Harvard Lyman- α
total water

ALIAS multipass laser IR
total water isotopes







Instruments

AVE-WIIF

- Lyman α water vapor
- **Lyman α total water**
- JPL tunable diode laser hygrometer (JLH)
- **Water vapor isotopes photodissociation laser-induced fluorescence (Hoxotope)**
- Water vapor isotopes integrated cavity output spectrometer (ICOS)
- ALIAS water isotopes

CRAVE

- Lyman α water vapor
- JPL tunable diode laser hygrometer (JLH)
- Water vapor isotopes integrated cavity output spectrometer (ICOS)
- ALIAS water isotopes
- **NOAA frostpoint (WB-57)**
- **Cryogenic frostpoint hygrometer (CFH)**
- **Balloon laser hygrometer**



What do these and other intercomparison plots tell us about instrument performance?

$$H_2O_{meas} = C(Temperature, Pressure, H_2O_{amb}) \times$$

$$\text{Signal}(Temperature, Pressure, H_2O_{amb})$$

$$H_2O_{meas} = C(Temperature, Pressure, H_2O_{amb} \pm H_2O_{inst} \pm H_2O_{???}) \times$$


$$\text{Signal}(Temperature, Pressure, (H_2O_{amb} \pm H_2O_{inst} \pm H_2O_{???}))$$

C is the instrument calibration factor.


H_2O_{amb} is the true water vapor mixing ratio

H_2O_{inst} is the sum of contaminant water vapor from either the instrument or aircraft.


Water vapor during AVE-WIIF




Harvard Lyman alpha and ICOS intercomparisons



Harvard Lyman alpha and ICOS intercomparison






Goal: To compare WB57 water vapor with CFH during the entire CRAVE mission.

Approach: Adjust JLH to agree with Lyman α during in situ part of mission (January 30 – February 11), then use JLH* during the remote part of the mission (January 14 – January 27) as a surrogate for Lyman α ,


Where $JLH^* = 0.755 \cdot JLH + 1.1$

Comparison: Lyman α and ICOS with CFH during CRAVE in situ and JLH* and ICOS for CRAVE remote.


Intercomparison of in situ and CFH during CRAVE in situ: 20060201, 20060209, 20060211




Intercomparison of in situ and CFH during CRAVE remote: 20060114,20060117,20060119,20060121,20060125




Intercomparison of in situ and CFH during CRAVE




Intercomparison of in situ and CFH during CRAVE



Intercomparison of in situ and CFH during CRAVE




Intercomparison of in situ and CFH during CRAVE




Conclusions

- As in AVE-WIIF, the overall agreement between Harvard water vapor instruments during CRAVE was very good.
- Comparisons between in situ water vapor on the WB57 and the CFH instrument illustrate systematic differences that increase significantly at low water vapor.
- Missions that provide the opportunity for careful water intercomparisons continue to be very useful and need to continue.
- Laboratory intercomparisons with low water vapor mixing ratios need to be carried out to help determine the source of this discrepancy.


Balloon TDL CRAVE sonde data




CFH CRAVE sonde data







Balloon TDL CRAVE sonde data



CFH CRAVE sonde data



AVE-WIIF, 20050703, 20050705, 20050707



The salient points of the operation, calibration, and in-flight validation of the Harvard water vapor are:

- Calibrations are carried out at a range of pressures and water vapor mixing ratios that are traceable to both the vapor pressure of water over liquid at room temperature and the absorption cross section of water vapor at Lyman- α (121.6 nm)
- In-flight measurements are carried out over a range of flow velocities (typically 20–80 m/s) to validate insensitivity to wall effects.
- ΔH_2O , measured by fluorescence using the laboratory calibration is cross-checked against ΔH_2O measured by dual path absorption. This cross-check validates the applicability of the laboratory calibration to in-flight conditions.
- During recent campaigns, agreement between the Harvard water vapor and total water instruments (in clear air) validates the insensitivity of the water vapor laboratory calibration to the in-flight temperature of the detection axis.